


XT Purification Kit Product Focus

Successful cycle sequencing reactions would not bring high-quality sequencing results without an effective sequencing reaction purification method. Although the traditional ethanol precipitation is cost-effective, the multiple step purification procedure needs extra caution from operators to prevent sample loss.

XT Purification Kit is a high-performance cost-effective alternative to XTerminator Purification kit. The undesired components after sequencing reactions such as unused dyes and dNTPs are removed by the purification resin without washing steps, reducing the sample loss possibilities. While saving is important for the service providers, not losing customers' samples is also essential for delivering high-quality results.

In addition to XT Purification Kit, we also provide more cost-effective Sequencing Reaction Cleaning Beads. For this product, we optimized the wash steps to minimize sample loss.

XT Purification Kit has different packaging sizes to meet the needs for different users.



XT Purification Kit

Catalog Numbers	Unit Size
8500050	50preps
8502010	1000preps
8502050	5000preps

Sequencing Reaction Cleaning Beads

Catalog Numbers	Unit Size
8005070	7ml
8006028	28ml
8006112	112ml
8007050	500ml